













Sterile Battery Operated Disposable Cautery Pens

These devices are able to achieve rapid hemostasis by precise cauterisation of the small blood vessels in a convenient and low cost package. They are extremely easy to handle being powered by an internal long life battery.

The range consists of five models which differ from each other in their working temperatures (1200°C, 800°C and 600°C) and the shape, thickness and length of the cautery tip. The range allows a wide variety of small incisions, resections or recesses to be dealt with effectively.

product no.	F7244	F7266	F7255	F7288	F7234	F7277
						
tip length	15mm	15mm	15mm	15mm	75mm	125mm
tip configuration						
temperature	1200	1200	600	800	1200	1200
power supply	3.0 V	3.0 V	1.5 V	3.0 V	3.0 V	3.0 V
weight	59 g	59 g	45 g	59g	72 g	78 g
unit price	£6.50	£6.50	£6.50	£6.92	£8.03	£8.03
pack price (10)	£61.70	£61.70	£61.70	£65.80	£76.30	£76.30

Application Guide

CARDIOVASCULAR & THORACIC SURGERY	Primary Choice	Secondary Choice
valve replacement	7234	7244
carotid endarterectomy	7244	7234
femoral popliteal artery bypass	7234	7244
patent ductus arteriosus	7234	7244
fistula formation	7244	7234
pacemaker implants	7244	7234
arteriovenous shunt	7234	7244
cardiac catheterism	7244	7234
PROCTOLOGY	Primary Choice	Secondary Choice
surgery on rectal condyloma	7277	7266
hemorrhoidectomy	7277	7244
treating of perineal laceration	7277	7244
GENERAL SURGERY	Primary Choice	Secondary Choice
colon resection	7277	7234
GYNAECOLOGY	Primary Choice	Secondary Choice
salpingectomy	7244	7234
abdominal hysterectomy	7244	7234
mastectomy	7244	7234
vulvectomy	7244	7234
vaginoplasty	7244	7234
GYNECOLOGY-Day Surgery	Primary Choice	Secondary Choice
incision and drainage of abscess	7277	7234
removal of vulvar condyloma acuminata	7277	7234
cervicitis polypectomy	7277	7234
treating of erosion and recurring cervicitis	7277	7234
excision of dysplasic tissue of little dimensions	7277	7234
NEUROSURGERY	Primary Choice	Secondary Choice
laminectomy	7277	7234
removal of tumors	7266	7244
acoustic neuroma	7244	
ORTHOPEDIC	Primary Choice	Secondary Choice
arthroplasty	7266	7234
knee surgery	7266	7234
hygroma removal	7234	
bone graft procedures	7234	
OTOLARYNGOLOGY	Primary Choice	Secondary Choice
thyroidectomy	7234	7244
control of eplstaxis	7277	7234
Polypectomy	7234	7277
tumor removal	7244	
carbuncle surgery	7266	7244
tracheostomy	7244	
excision of parotid gland	7266	7244
rhinoplasty	7244	
sinus surgery	7244	
nasal septal surgery	7244	
OBSTETRICS	Primary Choice	Secondary Choice
episiotomy repair	7244	
cesarean section	7244	
UROLOGY	Primary Choice	Secondary Choice
hydrolectomy	7244	